

INTRODUCER SHEATH

SPECIAL SHEATH FLEX						
		Sheath size F	Length	GW	GW Length	GW tip
PHIS 101	Introducer sheath flex 5F 8 cm	5	8	0,021	45 cm	J
PHIS 103	Introducer sheath flex 5F 11 cm	5	11	0,021	45 cm	J
PHIS 105	Introducer sheath flex 6F 8 cm	6	8	0,021	45 cm	J
PHIS 107	Introducer sheath flex 6F 11 cm	6	11	0,021	45 cm	J
PHIS 109	Introducer sheath flex 7F 11 cm	7	11	0,021	45 cm	J
PHIS 111	Introducer sheath flex 8F 11 cm	8	11	0,021	45 cm	J
PHIS 113	Introducer sheath flex 9F 11 cm	9	11	0,021	45 cm	J
PHIS 115	Introducer sheath flex 10F 11 cm	10	11	0,021	45 cm	J
PHIS 121	Introducer sheath flex 5F 24 cm	5	24	0,035	45 cm	J
PHIS 123	Introducer sheath flex 5F 45 cm	5	45	0,035	145 cm	J
PHIS 125	Introducer sheath flex 5F 65 cm	5	65	0,035	145 cm	J
PHIS 127	Introducer sheath flex 5F 90 cm	5	90	0,035	145 cm	J
PHIS 131	Introducer sheath flex 6F 11 cm	6	11	0,035	45 cm	J
PHIS 133	Introducer sheath flex 6F 24 cm	6	24	0,035	45 cm	J
PHIS 135	Introducer sheath flex 6F 45 cm	6	45	0,035	145 cm	J
PHIS 137	Introducer sheath flex 6F 65 cm	6	65	0,035	145 cm	J
PHIS 139	Introducer sheath flex 6F 90 cm	6	90	0,035	145 cm	J
PHIS 141	Introducer sheath flex 7F 11 cm	7	11	0,035	45 cm	J
PHIS 142	Introducer sheath flex 7F 24 cm	7	24	0,035	45 cm	J
PHIS 143	Introducer sheath flex 7F 35 cm	7	35	0,035	145 cm	J
PHIS 144	Introducer sheath flex 7F 45 cm	7	45	0,035	145 cm	J
PHIS 145	Introducer sheath flex R 7F 45 cm	7	45	0,035	145 cm	J
PHIS 146	Introducer sheath flex 7F 65 cm	7	65	0,035	145 cm	J
PHIS 147	Introducer sheath flex 7F 80 cm	7	80	0,035	145 cm	J
PHIS 148	Introducer sheath flex C 7F 90 cm	7	90	0,035	145 cm	J
PHIS 149	Introducer sheath flex 7F 100 cm	7	100	0,035	145 cm	J
PHIS 150	Introducer sheath flex 8F 11 cm	8	11	0,035	45 cm	J
PHIS 151	Introducer sheath flex 8F 24 cm	8	24	0,035	45 cm	J
PHIS 162	Introducer sheath flex 8F 45 cm	8	45	0,035	145 cm	J
PHIS 163	Introducer sheath flex 8F 65 cm	8	65	0,035	145 cm	J
PHIS 164	Introducer sheath flex 8F 80 cm	8	80	0,035	145 cm	J
PHIS 165	Introducer sheath flex 8F 100 cm	8	100	0,035	145 cm	J
PHIS 171	Introducer sheath flex 9F 11 cm	9	11	0,035	45 cm	J
PHIS 173	Introducer sheath flex 9F 24 cm	9	24	0,035	45 cm	J
PHIS 175	Introducer sheath flex 9F 45 cm	9	45	0,035	145 cm	J
PHIS 177	Introducer sheath flex 9F 65 cm	9	65	0,035	145 cm	J
PHIS 178	Introducer sheath flex 9F 80 cm	9	80	0,035	145 cm	J
PHIS 179	Introducer sheath flex 9F 100 cm	9	100	0,035	145 cm	J
PHIS 180	Introducer sheath flex 10F 11 cm	10	11	0,035	45 cm	J
PHIS 181	Introducer sheath flex 10F 24 cm	10	24	0,035	45 cm	J
PHIS 182	Introducer sheath flex 10F 45 cm	10	45	0,035	145 cm	J
PHIS 183	Introducer sheath flex 10F 65 cm	10	65	0,035	145 cm	J
PHIS 184	Introducer sheath flex 10F 80 cm	10	80	0,035	145 cm	J
PHIS 185	Introducer sheath flex 10F 100 cm	10	100	0,035	145 cm	J
PHIS 191	Introducer sheath flex 11F 65 cm	11	65	0,035	145 cm	J
PHIS 192	Introducer sheath flex 11F 80 cm	11	80	0,035	145 cm	J

SHEATH AORTIC						
		Sheath size F	Length	GW	GW Length	GW tip
PHIS 201	Introducer sheath AO 12F 28 cm	12	28	0,035	80	S
PHIS 203	Introducer sheath AO 14F 28 cm	14	28	0,035	80	S
PHIS 205	Introducer sheath AO 16F 28 cm	16	28	0,035	80	S
PHIS 207	Introducer sheath AO 18F 28 cm	18	28	0,035	80	S
PHIS 209	Introducer sheath AO 20F 28 cm	20	28	0,035	80	S
PHIS 211	Introducer sheath AO 22F 28 cm	22	28	0,035	80	S
PHIS 213	Introducer sheath AO 24F 28 cm	24	28	0,035	80	S
PHIS 215	Introducer sheath AO 26F 28 cm	26	28	0,035	80	S
PHIS 219	Introducer sheath AO 14F 64 cm	14	64	0,035	145	S
PHIS 221	Introducer sheath AO 16F 64 cm	16	64	0,035	145	S
PHIS 223	Introducer sheath AO 18F 64 cm	18	64	0,035	145	S
PHIS 225	Introducer sheath AO 20F 64 cm	20	64	0,035	145	S
PHIS 227	Introducer sheath AO 22F 64 cm	22	64	0,035	145	S
PHIS 231	Introducer sheath AO 26F 64 cm	26	64	0,035	145	S

SPECIAL SHEATH MICROPUNCTURE							
		Sheath size F	Length	GW	Needle		
PHIS 301	Introducer sheath micro 4F 8 cm	4	8	0,021	7 cm	45 cm	S
PHIS 303	Introducer sheath micro 5F 8 cm	5	8	0,021	7 cm	45 cm	S
PHIS 305	Introducer sheath micro 5F 11 cm	5	11	0,021	7 cm	45 cm	S
PHIS 307	Introducer sheath micro 5F 24 cm	5	24	0,035	7 cm	60 cm	S
PHIS 309	Introducer sheath micro 6F 8 cm	6	8	0,021	7 cm	45 cm	S
PHIS 311	Introducer sheath micro 6F 11 cm	6	11	0,021	7 cm	45 cm	S
PHIS 313	Introducer sheath micro 6F 24 cm	6	24	0,035	7 cm	60 cm	S
PHIS 315	Introducer sheath micro 7F 8 cm	7	8	0,021	7 cm	45 cm	S
PHIS 319	Introducer sheath micro 7F 24 cm	7	24	0,035	7 cm	60 cm	S
PHIS 321	Introducer sheath micro 8F 11 cm	8	11	0,021	7 cm	45 cm	S
PHIS 323	Introducer sheath micro 8F 24 cm	8	24	0,035	7 cm	45 cm	S
PHIS 325	Introducer sheath micro 9F 11 cm	9	11	0,021	7 cm	45 cm	S
PHIS 327	Introducer sheath micro 9F 24 cm	9	24	0,035	7 cm	45 cm	S

RADIAL SHEATH						
			Length	GW	GW length	GW TIP
PHIS 401	Introducer sheath radial 4F 7 cm	4	7	0,018	45cm	S
PHIS 403	Introducer sheath radial 5F 7 cm	5	7	0,018	45cm	S
PHIS 405	Introducer sheath radial 6F 7 cm	6	7	0,018	45cm	S
PHIS 407	Introducer sheath radial 7F 7 cm	7	7	0,018	45cm	S
PHIS 409	Introducer sheath radial 4F 11 cm	4	11	0,018	45cm	S
PHIS 411	Introducer sheath radial 5F 11 cm	5	11	0,018	45cm	S
PHIS 413	Introducer sheath radial 6F 11 cm	6	11	0,018	45cm	S
PHIS 415	Introducer sheath radial 7F 11 cm	7	11	0,018	45cm	S
PHIS 419	Introducer sheath radial 5F 16 cm	5	16	0,018	45cm	S
PHIS 421	Introducer sheath radial 6F 16 cm	6	16	0,018	45cm	S
PHIS 423	Introducer sheath radial 7F 16 cm	7	16	0,018	45cm	S
PHIS 425	Introducer sheath radial 4F 23 cm	4	23	0,018	45cm	S
PHIS 427	Introducer sheath radial 5F 23 cm	5	23	0,018	45cm	S
PHIS 431	Introducer sheath radial 7F 23 cm	7	23	0,018	45cm	S

FEMORAL SHEATH						
		Size	Length	GW	GW Legnth	GW TIP
PHIS 501	Introducer sheath femoral 4F 11 cm	4	11	0,035	50 cm	S
PHIS 503	Introducer sheath femoral 5F 11 cm	5	11	0,035	50cm	S
PHIS 505	Introducer sheath femoral 6F 11 cm	6	11	0,035	50 cm	S
PHIS 507	Introducer sheath femoral 7F 11 cm	7	11	0,035	50 cm	S
PHIS 509	Introducer sheath femoral 8F 11 cm	8	11	0,035	50cm	S
PHIS 511	Introducer sheath femoral 4F 23 cm	4	23	0,035	50 cm	S
PHIS 513	Introducer sheath femoral 5F 23 cm	5	23	0,035	50 cm	S
PHIS 515	Introducer sheath femoral 6F 23 cm	6	23	0,035	50cm	S

FEMORAL SHEATH PEEL AWAY						
		Size	Length	GW	GW Legnth	GW TIP
PHIS 601	Introducer sheath peel away 4F 11 cm	6	13	0,038	60 cm	S
PHIS 603	Introducer sheath peel away 5F 11 cm	7	13	0,038	60 cm	S
PHIS 605	Introducer sheath peel away 6F 11 cm	8	13	0,038	60 cm	S
PHIS 607	Introducer sheath peel away 7F 11 cm	9	13	0,038	60 cm	S
PHIS 609	Introducer sheath peel away 8F 11 cm	10	13	0,038	60 cm	S
PHIS 611	Introducer sheath peel away 4F 23 cm	10,5	13	0,038	60 cm	S
PHIS 613	Introducer sheath peel away 5F 23 cm	11	13	0,038	60 cm	S
PHIS 615	Introducer sheath peel away 6F 23 cm	12	13	0,038	60 cm	S
PHIS 619	Introducer sheath peel away 6F 23 cm	6	23	0,038	60 cm	S
PHIS 621	Introducer sheath peel away 7F 23 cm	7	23	0,038	60 cm	S
PHIS 623	Introducer sheath peel away 8F 23 cm	8	23	0,038	60 cm	S
PHIS 625	Introducer sheath peel away 9F 23 cm	9	23	0,038	60 cm	S
PHIS 629	Introducer sheath peel away 10,5F 23 cm	10,5	23	0,038	60 cm	S
PHIS 631	Introducer sheath peel away 11F 23 cm	11	23	0,038	60 cm	S